

Adosy_ste_setup_BBFO Topspin 3.6 600MHz. Fu Chen October 2021

Usage: Optimize critical parameters for DOSY experiment for solution without solvent suppression at BBFO probe.

=====

Parameters to be optimized: P30, d20, RG and NS.

For standard D2O/Cu test sample, sine.100 gradient used, P30 1.0ms, d20 0.06s; D1 7s; NS 1

*****WARNING:

DONOT increase P30 more than 2.5ms.

=====

General Instruction:

1. Create experiment #1: Set stable temperature (type edte), tune and shim well.
 - Run 1H NMR (Ah1): obtain o1p and sw
 - Remark: For 100% water, detune the probe by 10MHz, and calibrate 90.
2. Calibrate proton 90 degree pulse at #1. Obtain P1 and PL1 value for your sample.
 - Automation: type pulsecal
 - Manually: RPAR Ah1_90_water
3. Create experiment #2:
 - rpar Adosy_ste_setup_BBFO
 - getprosol
4. For Staff only:
 - Check the record for gradient calibration (gradpar cf with edhead gradient history).
5. Update o1p, sw and P1 (from step 1-2), NS 1, D1 7s to #2
 - Set d20 (default 60ms, typical 50-150ms)
 - Check to use the following default values as a starter {gpz6 2; gpz7 -17.13 (fixed); gradient pulse shape: SINE.100}
6. Optimize: P30 from 1ms (typical 1-2ms sample dependent)
 - Set a value of p30 between 1.3-1.5ms for water,
 - Demo sample (H2O/D2O with Cu+2). uses 1.0ms.
 - rga. and zg to collect and phase the first spectrum properly.
7. Do three measurements vary parameter gpz6:

- At #2, type i or iexpno to create #3
 - Change gpz6 to 95 at #3
 - At #2, type multizg, then type 2
 - At #2 and #3, type efp;apk
 - Compare peak intensities for #2 and #3 using multi-display; signals should decrease approx. 50 ratio.
8. Adjust P30 values:
- If Signal decay too rapidly, return to #2 and #3, decrease the P30 value by 0.1ms; repeat multizg at step 7.
 - If Signal decay too slowly, return to #2 and #3, increase the P30 value by 0.1ms; repeat multizg at step 7.
 - Remark: P30 should not exceed 2.5ms. If P30 is out of its range, go to the next step.
9. If step 8 fails, adjust d20 by 0.1ms and repeat steps 7 & 8, then step 9 as needed.
- d20 (default 60ms, typical 50-150ms)
 - The larger the P30 or/and D20, the faster the decay in intensities.
10. Optimize NS for your sample:
- Set GPZ6 set to 95% with NS = 8
 - zg to collect a spectrum; check the sufficient intensities of the signal of interests by increasing NS = 8